NAME:	PHONE:	
HOSPITAL/ORGANIZATION:	E-MAIL:	
DEALER/REP:	SALES REP:	QUOTE:



WWW.PHSMEDICALSOLUTIONS.COM 800-743-7738 • FAX: 605-882-8398

## **E-EPD** *Electric Epidural Positioning Device*

## **STANDARD FEATURES:**

- Load Capacity: For use with patients up to 600 lbs.
- 28"W x 29"D x 50-60"H
- White powder coat base
- 180° Adjustable face rest
- Prepositional pillow for use during cervical, lumbar and thoracic placement
- Adjustable arm rest
- Arm rest size: 21"W x 15.5"D
- Arm rest adjustment: 6"
- Adjustable torso support
- Torso cushion size: 6.25"W x 5.75"H
- Torso cushion adjustment: 6"
- Single hand control adjusts both midsection and foot rest
- (2) Rechargeable batteries last up to 100 cycles each
- Adjustable foot rest
- Medical grade profile locking casters
- Black upholstery
- Includes 1 Box (50) Sani-Cover® Fitted Disposable Face Rest Covers



Based on a design by an anesthesiologist, the E-EPD positions patients correctly and comfortably by encouraging cervical, thoracic and lumbar flexion while maintaining a solid and stable position. A new curved design, for use with patients up to 600 lbs., accommodates patients of a variety of body sizes and weights. Electronic controls adjust the E-EPD to the correct height and foot position for each patient, while the easy mechanical lock and release system makes chest, arm and head adjustments quick and effortless.

The E-EPD also offers these other benefits:

- It's intuitive to use and easy to set up, clean and move between rooms.
- Convenient, disposable face rest covers make for fast and sanitary changeovers between patients.
- A single foot pedal easily locks or unlocks all 4 casters quickly and effortlessly.
- The battery charges fully in 6 hours and lasts between 40 to 100 cycles (based on load).

## STANDARD E-EPD:

☐ E-EPD

Electric Epidural Positioning Device

Qty \_\_\_\_

## **OPTIONS:**

☐ **FRC-50-4221 -** Sani-Cover® Fitted Disposable Face Rest Covers (50 Covers per Box)

Qty \_\_\_\_